


Spinsolve 60 Carbon


60 MHz carbon benchtop NMR

The Spinsolve 60 Carbon benchtop NMR spectrometer gives you a remarkable 60 MHz NMR spectrometer in a compact benchtop instrument. The spectrometer is capable of a wide range of ¹H, ¹⁹F, and ¹³C NMR experiments. The software is beautiful and easy for anyone to use. You can also quickly automate running multiple experiments with flexible scripting.







Spectra from the Spinsolve 60 Carbon benchtop NMR




2 M Ibuprofen, 1D proton, single scan, 10 seconds



2 M Ibuprofen, HSQC-ME, ~1 hour




Neat Diethyl phthalate. 1D Carbon plus three DEPT spectra, total time 10 minutes



2 M Ibuprofen, HMBC, ~2 hours

Software

The Spinsolve software is beautifully simple and easy to use, with a clean and intuitive user interface.



Simple menu structure


- 1. Click to choose nucleus
- 2. Click to choose experiment
- 3. Click Start (watch status on the progress bar)
- 4. Click any processing you wish to apply



Automate experiments with scripts

A range of script templates are provided for easy user modification. Scripts are displayed graphically to provide a clear picture of the sequence.

This script runs a sequence of a 1D proton, a 1D carbon, a DEPT, an HSQC-ME, and an HMBC.







Specifications

• Nuclei: 1H, 19F, 13C

• Operating frequency: 60 MHz (1H)

• 1H 50% Linewidth: < 0.5 Hz

• 1H 0.55% Linewidth: < 20 Hz

• ¹H Sensitivity: >120:1 for 1% Ethyl Benzene

• Operating Temperature Range: 20° C to 25° C (68° F to 77° F)

• Dimensions: 58 x 43 x 40 cm (23" x 17" x16")

• Weight: 60 kg (120 lb)

Stray Field: < 2 G all around systemPower requirement: 110-240V AC



Pulse sequences available on the Spinsolve 60 Carbon spectrometer

Proton	Fluorine	Carbon
1D	1D	1D
Paramagnetic	Paramagnetic	DEPT
2D COSY	2D F - COSY	HETCOR
2D TOCSY	2D F - JRES	HMBC
2D JRES	2D FH - COSY	HMQC
T_1,T_2	Reaction Monitoring	HSQC
Reaction Monitoring		HSQC-ME

Other sequences may be available, contact Magritek for details.

Contact us now for a quote, to request a demo or to measure your samples

Email: sales@magritek.com

Website: www.magritek.com/contact-us

GERMANY +49 241 9278 7270 **UNITED STATES** +1 855 667 6835 **UNITED KINGDOM** +44 7468 529 615 **NEW ZEALAND** +64 4 477 7096

For a complete listing of our global offices and distributors visit: www.magritek.com/about-us/distributors/